

ADA Lap System

ADA Lap (Pronounced add-a-lap) Solutions provide a universal support for wheelchair users.



The ADA Lap™ simply slides under your power or manual wheelchair cushion to provide a versatile base for attachments.

ADA 01 - Base Unit

The easy-to-use base simply slides under the wheelchair cushion.



It is held in place by the weight of the user. Height adjustment is by a quick release clamp and the vertical pole is easily removable for transfer.

Application

The ADA Lap Food Tray in use by Gerry Price



Before his workplace accident in 1995, Gerry worked in Engineering and Operations.

Since becoming a quadriplegic, he understands that to be alive and healthy, the disabled need to get out and be active in their communities.

Why it works:

- Easily removed for transfer
- Fits many wheelchairs, regardless of width
- Great for facilities (share between users)
- Unaffected wheel access for self propelling
- Does not wear out armpads
- Does not affect arm positioning
- Strong and lightweight aluminum
- Height and Angle adjustable

More Fine Products from Symmetric Designs

Headmaster Collar



TOT Collar



Wheelchair Headrests



ADA Lap™

Solutions for an Active Lifestyle



Universal Design to get you going!

With many attachment options, the ADA Lap system is a versatile wheelchair accessory mount that is easy to use, lightweight, low profile and compatible with any width of wheelchair.

The base slides under the wheelchair cushion and is held securely in place by the weight of the user.

MD The Axion Rotary Interface is a Medical Device

ADA Lap Accessories



The Lap Shopper

The Lap Shopper provides hands free convenience. This enables the shopper with limited dexterity to have both hands free to pick up merchandise from the shelf.

Fishing Rod Holders

Go fresh or salt water fishing with this great rod holder! The rod holder plugs into the extension tube of the ADA Lap Base.




Available in either fresh water or salt water models. For use with an Extension Pole. Compatible with Scotty Hardware.

ADA FW (Fresh Water Rod Holder)
ADA SW (Salt Water Rod Holder)



ADA Lap Accessories



Lexan Tray

This strong polycarbonate tray provides a firm transparent surface for writing.

ADA TR

Food Tray

With a lightly textured bottom and curved lip, this 11" x 14" sturdy tray will keep drinks and plates from sliding off. Also great for containing any small objects that may roll or slide on a flat surface.



ADA FT



SNAP Lap Desk

A small attache case that opens to contain your papers, pens, etc., yet doubles as a desk. Use with a non-slip surface for keeping control of slippery objects. No more juggling items on your lap!

ADA LD



Extension Poles

Available in either 12" or 18" lengths, extension poles have an insert for camera mounts/tripods, Umbrella holders, Angle Adapters and fishing rod holders.



ADA PL-12 ADA PL-18

ADA Lap Accessories

Angle Adapter

This versatile adapter comes with a removable top plate surfaced with velcro. It can be used to angle trays and other attachments. For use with an extension pole.



ADA AA

Camera Mount

The ADA Lap camera mount/tripod allows small height adjustments using a hand crank as well as allowing camera tilt function. This camera mount can be used on the ADA Lap or freestanding. For use with an extension pole when mounted.



ADA CM

Umbrella Connector

The ADA Lap Umbrella connector can be screwed into a pre-drilled hole in any wooden handled umbrella. Also available pre-installed with umbrella. For use with an extension pole.



ADA UC or UB